



BRASSELER
USA®

BY YOUR SIDE

ENDOSEQUENCE® BC SEALER™ & BC SEALER HIFLOW™


THE "STATE-OF-THE-ART" IN ENDODONTIC OBTURATION
IS NOW AVAILABLE FOR WARM TECHNIQUES!

Clinicians now have two distinct formulations of the same proven pre-mixed bioceramic sealer to choose from. Compared to standard BC Sealer, BC Sealer HiFlow exhibits a lower viscosity when heated and is more radiopaque, making it optimized for warm obturation techniques.

BC Sealer 2g

Order Number: 5017560U0

- Sealer-based technique
- Optimal for cold obturation



NEW!


BC Sealer HiFlow 1.5g

Order Number: 5026825U0

- Condensation-based technique
- Optimal for warm obturation
- Lower viscosity when heated
- Higher radiopacity



BIOCERAMIC BONDED OBTURATION

Standard gutta percha can be used with BC Sealer, but for a true, gap-free seal, BC Points are recommended. BC Points are impregnated and coated with bioceramic particles to allow for bonding with BC Sealer. The benefit of using BC Sealer and BC Points is that three-dimensional bonded obturation can be achieved at room temperature. However, if a warm vertical technique is preferred it is recommended to use BC Sealer HiFlow with BC Points 150 Series along with BC Pellets at 150°C - 220°C.



BC SEALER SETTING/BONDING REACTION

BC Sealer is a calcium silicate, calcium phosphate based sealer. The calcium silicate combines with water to form calcium hydroxide (high pH). The calcium hydroxide then combines with calcium phosphate to form hydroxyapatite.



COLD BONDED OBTURATION

- BC Sealer
- BC Points



EndoSequence BC Points™
Individual Packs (60/Pack)

✓ *Stiffer to drive sealer*

Tip ID	Tip Size	.04 Taper	.06 Taper
○ White	15	5019508U0	5019520U0
● Yellow	20	5019509U0	5019521U0
● Red	25	5019510U0	5019522U0
● Blue	30	5019511U0	5019523U0
● Green	35	5019512U0	5019524U0
● Black	40	5019513U0	5019525U0
○ White	45	5019514U0	5019526U0
● Yellow	50	5019515U0	5019527U0
● Red	55	5019516U0	-
● Blue	60	5019517U0	-
● Green	70	5019518U0	-
● Black	80	5019519U0	-

WARM BONDED OBTURATION

- BC Sealer HiFlow
- BC Points 150 Series & BC Pellets



EndoSequence BC Points™ 150 Series
Individual Packs (60/Pack)

✓ *Melt at 150°C*

Tip ID	Tip Size	.04 Taper
○ White	15	5024583U0
● Yellow	20	5024584U0
● Red	25	5024585U0
● Blue	30	5024586U0
● Green	35	5024587U0
● Black	40	5024588U0
○ White	45	5024589U0
● Yellow	50	5024590U0

BC Pellets (100/Pack)

✓ *Melt at 150°C*

Order No. 5023074U0



Visit our website at BrasselerUSA.com

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

