




# KontrolFlex™ Rotary Files

## CLINICAL TIPS

- Always reach full working length to at least a #15/.02 prior to use of rotary files.
- 500-600 RPMs.
- 1.8-2.3 Ncm.
- Use a feather-light touch taking the file to engagement (not resistance).
- Use the Single Stroke and Clean Technique. Clean the flutes after each gentle engagement.
- Irrigate between each rotary file.
- Obturate the canals using BC Sealer™ and the corresponding BC Point™.

### KONTROLFLEX™ CROWN DOWN TECHNIQUE

- For crown down, move to the next smallest file after 4-5 strokes of single engagement.



Patent Numbers: 8911573, 9005377

Order online at [BrasselerUSA.com](http://BrasselerUSA.com) or call 800.841.4522.  
In Canada call 800.363.3838.



**BRASSELER**  
USA® BY YOUR SIDE™




# KontrolFlex™ Rotary Files

## CLINICAL TIPS

- Always reach full working length to at least a #15/.02 prior to use of rotary files.
- 500-600 RPMs.
- 1.8-2.3 Ncm.
- Use a feather-light touch taking the file to engagement (not resistance).
- Use the Single Stroke and Clean Technique. Clean the flutes after each gentle engagement.
- Irrigate between each rotary file.
- Obturate the canals using BC Sealer™ and the corresponding BC Point™.

## KONTROLFLEX™ VARIABLE TAPER TECHNIQUE

- Begin sequence with #19/.07 taper taking it to only 2-3 mm into the canal. Do not take to working length.
- Proceed using variable taper technique using the files in a progressive series. Move to the next file after the current file reaches working length (see chart).



Patent Numbers: 8911573, 9005377

Order online at [BrasselerUSA.com](http://BrasselerUSA.com) or call 800.841.4522.  
In Canada call 800.363.3838.



**BRASSELER**  
USA® BY YOUR SIDE™