



BRASSELER
USA®

ENDOSEQUENCE® CM CONTROLLED MEMORY FILES

*The Perfect Balance of
Strength and Flexibility*

BY YOUR SIDE




ENDOSEQUENCE® CM

THE PATENTED ENDOSEQUENCE® DESIGN

The original EndoSequence file was introduced over 15 years ago, and despite the introduction of countless competitive files, EndoSequence remains one of the most popular files available today. The patented design and quality manufacturing set the EndoSequence file apart and make it one of the most efficient cutting files on the market.

WHAT YOU LOVE ABOUT ENDOSEQUENCE®, MADE BETTER

EndoSequence CM Files are manufactured with an enhanced nitinol metallurgy and heat treatment process that allows for cutting edge fidelity and superior resistance to cyclic fatigue. EndoSequence CM provides clinicians with a more flexible option and an instrument that is significantly more resistant to separation.



ENDOSEQUENCE® CM PRODUCT OFFERING

(4/PACK)

EndoSequence CM Files are provided in sterile blister packs with four instruments per pack.

STERILE R



File/ Tip ID	Tip Size	.04 Taper		.06 Taper	
		21 mm	25 mm	21 mm	25 mm
△ White	15	5026609U0	5026617U0	5026633U0	5026639U0
▲ Yellow	20	5026610U0	5026618U0	5026634U0	5026640U0
▲ Red	25	5026611U0	5026619U0	5026635U0	5026641U0
▲ Blue	30	5026612U0	5026620U0	5026636U0	5026642U0
▲ Green	35	5026613U0	5026621U0	5026637U0	5026643U0
▲ Black	40	5026614U0	5026622U0	5026638U0	5026644U0
△ White	45	5026615U0	5026623U0	—	—
▲ Yellow	50	5026616U0	5026624U0	—	—

File/ Tip ID	Tip Size	Variable Taper			
		19 mm	21 mm	25 mm	31 mm
▲ N/A	SX/19	5027891U0	—	—	—
▲ Purple	S1/18	—	5027892U0	5027899U0	5027906U0
△ White	S2/20	—	5027893U0	5027900U0	5027907U0
▲ Yellow	F1/20	—	5027894U0	5027901U0	5027908U0
▲ Red	F2/25	—	5027895U0	5027902U0	5027909U0
▲ Blue	F3/30	—	5027896U0	5027903U0	5027910U0
▲ Black	F4/40	—	5027897U0	5027904U0	5027911U0

PREMIUM BRASSELER QUALITY


EndoSequence CM files are ultra-premium NiTi files that adhere to extremely tight quality control specifications. Each step of the manufacturing process is meticulously controlled including our exclusive electropolishing process which provides the cleanest and sharpest cutting edge available on the market. EndoSequence CM files are provided in sterile blister packs to ensure that you are using the cleanest and most efficient cutting file for every procedure.

ELECTROPOLISHED SURFACE

- Eliminates surface imperfections that weaken other files
- Increases file sharpness




EndoSequence® CM File




Vortex Blue™ File

PRECISION TIP

- Non-active at the tip, fully active precisely at 1 mm



EndoSequence® CM File



Vortex Blue™ File


ENHANCED NITINOL METALLURGY

- Superior resistance to cyclic fatigue
- Allows for pre-bending for difficult access
- Deformation occurs prior to separation
- Cutting edge fidelity maximized at body temp



PATENTED ACP (ALTERNATING CONTACT POINT™) DESIGN

- Non-landed, reduced-drag design
- Centers file within canal without radial lands
- Prevents self-threading "screw in" to canal
- Enhances debris removal
- Maximizes cutting efficiency



NO RADIAL LANDS

- Flexibility is maximized
- Torque is minimized
- Reduced cross section thickness
- Sharp cutting edge



3 point engagement



2 point engagement



1 point engagement




No engagement

STERILE PACKAGING


- Eliminates cross-contamination
- Improves practice efficiency

STERILE R


Vortex Blue™ is a registered trademark of Dentsply Sirona, Inc.



DIAMONDS



CARBIDES



POLISHERS



PROCEDURE SYSTEMS



BIOCERAMICS



Since 1976, our absolute focus has been to develop products and provide services that support the practice of core operatory procedures. Today, Brasseler USA® is the premier dental instrumentation company in North America offering the most extensive selection of dental and surgical instrumentation under one brand.



HANDPIECES



ENDODONTICS



SMALL EQUIPMENT



HYGIENE



HAND INSTRUMENTS



ORAL SURGERY



LABORATORY INSTRUMENTS



Visit our website at BrasselerUSA.com

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

