EndoSequence® BC Liner™ and BC Liner™ Blue Instructions for Use

INDICATIONS FOR BC LINER™

Recommended as a bioactive base and liner for Class I, II, III and V restorations where there is no pulpal involvement, for indirect pulp capping, and for use with all composite and amalgam restorations.

CONTRAINDICATIONS FOR BC LINER™

Not indicated for direct placement on the exposed pulp. See instructions for use.

BC LINER™ INSTRUCTIONS FOR USE

1. Set curing light to 20-second low intensity setting.
2. Isolate and prepare tooth to receive a restoration. Etching and bonding agents are not required.
3. Place appropriate pulp protection, if indicated.
4. Lightly dry removing excess moisture with high volume evacuation, compressed air or cotton pellet. Do not desiccate tooth.
5. Place a mix tip on the BC Liner™ syringe. To ensure an even mix of base and catalyst, dispense 1-2 mm of material onto a mixing pad and discard this material.
6. Open Sandwich Technique: Apply BC Liner™ to prepared surfaces and extend to the enamel-cavo surface margin. Massage into the dentin for 20 seconds and light cure using 20-second low intensity light setting. Initial self cure setting time is 5 minutes. Continue with step 8.
7. Closed Sandwich Technique: Apply BC Liner™ to prepared surfaces. Do not extend the material over the enamel margins. Massage into the dentin for 20 seconds and light cure using 20-second low intensity light setting. Initial self-cure setting time is 5 minutes.
8. Complete the restoration with your preferred adhesive and composite bonding system. Finish and polish as usual.

Note: Allowing dual cure material to self-cure for 20-30 seconds before light curing mitigates polymerization stress and exothermic reaction. Always use 20-second low intensity light setting. Do not exceed manufacturer's recommended curing time. Additional curing may generate excessive heat damaging to the pulp.

Note: Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient. For additional information, refer to: <http://www.fda.gov/Medical-Devices/dental-devices/multiple-use-dental-dispenser-devices>

HOW TO USE THE AUTOMIX SYRINGE

1. Remove cap. If necessary bleed the syringe so that base and catalyst are at the orifice of the syringe barrels. Place a mixing tip on the automix syringe.
2. To ensure an even mix of base and catalyst, dispense 1-2 mm onto a pad and discard this material.
3. Dispense material directly onto the tooth or into the restoration.
4. Discard mixing tip. Recap syringe. Do not cross-contaminate base and catalyst.

APPLY BC LINER™ TO A DRY TOOTH SURFACE, BUT DO NOT DESICCATE THE TOOTH

Using high volume evacuation, compressed air and/or a cotton pellet, dry and remove all external moisture from the prepared tooth surface. Do not desiccate the tooth, which naturally contains a small amount of water. Ceramic, metal, resin desensitized, prehybridized, and cured composite surfaces should be dry.

CAUTION

Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves. Do not exceed manufacturer's recommended curing time for the light you are using.

EndoSequence® BC Liner™ and BC Liner™ Blue Instructions for Use

STORAGE AND HANDLING

- Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.
- Shelf life of unopened product: 2 years from date of manufacture.
- Re-cap immediately after use.

REPORTING OF SERIOUS INCIDENTS

Any serious incidents involving the device should be reported to the manufacturer and the Competent Authority of the Member State in which the incident occurred, or other appropriate regulatory body. Reports of serious incidents should be sent to: Brasseler U.S.A. Dental, LLC One Brasseler Boulevard Savannah, Georgia 31419 USA Telephone: 1-800-841-4522.

Note: Refer to the Safety Data Sheet for a list of materials present in this product.

CURING PROPERTIES	BASE/LINER
Light cure setting time:	20 seconds <i>(See note in instructions)</i>
Working time at room temperature:	90 seconds
Initial self-cure setting time (mouth temperature):	≤ 5 minutes

Glossary of Symbols: www.Brasselerusadental.com/resources

Manufactured for
Brasseler U.S.A. Dental, LLC
One Brasseler Blvd.
Savannah, GA 31419 USA
Telephone 800-841-4522

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Brasseler U.S.A. Dental, LLC One Brasseler Boulevard Savannah, Georgia 31419 États-Unis Téléphone : 1-800-841-4522.