Brasseler Core Build-Up Material Dual-Cure with Fluoride

BRASSELER USA®

Store Under Refrigeration

Brasseler Core Build-up Material is a creamy, radiopaque core build-up and posterior restorative material with excellent syringeability and control. Brasseler Core Build-up Material is dual-cure and is available in White with Fluoride and Enamel with Fluoride.

INTENDED USE: Composite core build-up material. Use for the cementation of posts and pins, as a composite for crown build-ups, repairing and preventing split roots and for integral post and core procedures in one visit. Eliminates the need for a cast post. Enamel shade may also be used as a posterior restorative material.

CAUTION: Wear protective gloves while using this product.
CAUTION: Wear eye protection while using this product.
CAUTION: Brasseler Core Build-up Material cannot be used in children, pregnant or breast-feeding women.

DIRECTIONS FOR BRASSELER CORE BUILD-UP MATERIAL

Important: Bring Core Build-up Material to room temperature 30 minutes before using to ensure proper set time.

1. Clean and isolate the teeth.
2. If you are using a post, it is important that you etch the tooth/root surface after you size and fit the post per etchant's manufacturer's instructions.
3. Apply dental adhesive of your choice per manufacturer's instructions,however, ensure final curing of Brasseler Core Build-up Material cement by curing light.

Note: If you are using a light-cure bonding agent, it is important that you apply a coupling agent after the application of the bonding agent to ensure bonding compatibility. Then, gently air-dry the surface.

A coupling agent is recommended for use with self/dual-cure composites when using any light cure bonding agent that you are unable to reach with a curing light in order to activate.

4. Auto Mixing:
 a. Align the straight edge of the auto-mixing housing with the syringe flange.
 b. Push the mixing tip onto the syringe and turn 90° clockwise until it stops.
 c. Firmly attach the intraoral tip to the mixing tip.

Important: You must extrude a pea-sized amount of material after placing the mixing tip on the syringe - discard this material. This is critical to ensure that Brasseler Core Build-up Material sets properly. Repeat this each time a mixing tip is placed on the syringe.

5. Position the intraoral tip directly into the canal/preparation and dispense Brasseler Core Build-up Material.
6. After Brasseler Core Build-up Material has cured, prepare the tooth for final restoration.

Additional auto-mixing and intraoral tips are available.
Contact a Brasseler Sales Representative for ordering information.

DIRECTIONS FOR TOOTH-COLORED POSTERIOR RESTORATIVES

Brasseler Core Build-up Material Enamel Shade is also indicated for use as a posterior restorative. Brasseler Core Build-up Material has the ideal compressive strength and wear resistance for long-lasting posterior restorations. The syringeable formula and delivery system make placement faster and easier than conventional posterior composites.

1. Select the Enamel Shade of Brasseler Core Build-up Material.
2. Follow Steps 1–6, DIRECTIONS FOR CORE BUILD-UP.
3. Using a small 12-fluted football fine diamond bur, create anatomy details.
4. Finish and glaze: Use a fine diamond and a 12/30 fluted bur to finish the composite. Rinse and air-dry the tooth composite surface. Check occlusion and modify as needed. Apply a light-cured clear resin to add natural luster and additional wear resistance. Light-cure clear resin.

Due to variations in the performance characteristics of light curing units ALWAYS bench test restorative materials before use in vivo. Curing test rings are provided for this purpose.

a. Fill the 2mm deep well of the test ring and level material.
b. Position the light transmitting element perpendicular to and approximately 2mm–5mm above the top surface of the ring.

- With Sapphire PAC lights (all models), start with 5 second exposures.
- With Flashlite LED lights (all models), start with 10 second exposures. Approximately 2mm-5mm above the top surface of the ring.
- For all other curing lights – halogen, LED and others, refer to the manufacturer instructions. A minimum of 10–30 seconds is recommended.

c. Use a dental probe to scrape test the hardness of the top and bottom surfaces. The bottom surface should be as hard as the top surface.
d. If the bottom surface is not completely cured repeat steps (b) to (c). Repeat until the bottom surface is completely cured.
e. Maintain a log including material, shade and associated curing exposure time. Use the log to monitor system performance.

Note: If a cavity preparation is deep, curing exposure times must also be increased due to beam divergence and angular placement of the light transmitting element to the restoration. An incremental filling technique is recommended and each increment should be fully cured prior to applying additional layers.

General guidelines for curing light unit exposure times. See manufacturer's instructions. ALWAYS bench test restorative materials before use in vivo.

- Curing lights with power density greater than 800 mW/cm2 – Cure the buccal and lingual with 10 second exposures for each area.
- Curing lights with power density less than 800 mW/cm2 – Cure the buccal and lingual for 20 second exposures for each area.
- Curing lights with power density less than 300 mW/cm2 should not be used to cure.

STORAGE

Store under refrigeration.

Do not expose to temperatures exceeding 77° F (25° C).

Do not expose to direct sunlight.

Do not freeze.

SAFETY DATA SHEETS AVAILABLE AT BRASSELERUSADENTAL.COM

Manufactured for
Brasseler U.S.A. Dental, LLC
One Brasseler Blvd
Savannah, GA 31419 U.S.A.

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