




BRASSELER
USA®







BY YOUR SIDE

BRiO SHINE™



SINGLE-USE UNIVERSAL COMPOSITE POLISHING SYSTEM

- ▶ System adapts to all tooth surfaces
- ▶ Reduces time spent polishing
- ▶ Fewer instruments and steps

						
Order No.	BPTCO-1	BPTF-2	BCPCO-1	BCPF-2	FLBCO-1	FLBF-2
Head Size (1/10 mm)	040	040	060	060	100	100
Grit	Medium	Fine	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)						
Qty/Pack	30	30	30	30	30	30

Recommended speed: 4,000-6,000 rpm, Final Polish: 2,000-3,000 rpm

▶ Order online at BrasselerUSA.com or call **800.841.4522**.
In Canada call **800.363.3838**.



UNIVERSAL COMPOSITE
STEP 1

Pre-polishing composite restorations



UNIVERSAL POLISHING
STEP 2


Final, universal polisher that provides an excellent finish to all restorations



BRASSELER
USA®







BY YOUR SIDE

BRiO SHINE™



SINGLE-USE UNIVERSAL CERAMIC POLISHING SYSTEM

- ▶ System adapts to all tooth surfaces
- ▶ Reduces time spent polishing
- ▶ Fewer instruments and steps

						
Order No.	BPTCER-1	BPTF-2	BCPCER-1	BCPF-2	FLBCER-1	FLBF-2
Head Size (1/10 mm)	040	040	060	060	100	100
Grit	Medium	Fine	Medium	Fine	Medium	Fine
21 - RA - Right-Angle Latch (22 mm)						
Qty/Pack	30	30	30	30	30	30

Recommended speed: 4,000-6,000 rpm, Final Polish: 2,000-3,000 rpm

▶ Order online at BrasselerUSA.com or call 800.841.4522.
In Canada call 800.363.3838.



UNIVERSAL CERAMIC
STEP 1

Pre-polishing zirconia, lithium disilicate and porcelain/PFM restorations



UNIVERSAL POLISHING
STEP 2

Final, universal polisher that provides an excellent finish to all restorations