



BRASSELER
USA® SURGICAL

VERSADRIVER®

PNEUMATIC OSTEOTOME SYSTEM



Enhancing Surgical Precision & Efficiency in Hip & Knee Revisions

► KEY BENEFITS

- Increased control: Directional strikes, for precise separating control.
- Minimizes bone loss: Preserves anatomy during implant removal.
- Flexible blades: Conform to implant geometry for precise removal.
- Versatile: Press fit or cemented hip stems and knee tibial/femoral components.
- Reduced strain: Replaces chisels and mallets for less fatigue.
- Guided Control: Pistol grip design enables the surgeon to guide the blade by hand.





BRASSELER
USA® SURGICAL

VERSADRIVER®

PNEUMATIC OSTEOTOME SYSTEM



FEATURES

- Single or multi-strike separating actions
- Hand-piece & six Osteotome blades
- Up to 3900 strikes per minute
- Blade limited to 3mm stroke

PRODUCT CODE	ITEM DESCRIPTION
10-100	VersaDriver® Pneumatic Osteotome System
10-001	VersaDriver® Hand Piece
10-092	VersaDriver® Case
10-093	VersaDriver® Tray
10-095	Hose
10-069	Blade Extraction Slap Hammer
10-101	Flat Tip Blade 11mm x 50mm
10-102	Flat Tip Blade 15.5mm x 100mm
10-103	Round Tip Blade 11mm x 100mm
10-104	Round Tip Blade 15.5mm x 100mm
10-105	Flat Tip Blade 11mm x 100mm
10-112	Cement Removal Parallel V Tip Blade 15.5mm x 100mm



Call 800.535.6638 to order today!
Visit us online at BrasselerSurgical.com