

FORZA™ L50K

Uncompromising Power...

The Forza L50K Lab Micromotor delivers uncompromising power throughout its advanced design. The cutting edge technology incorporated in the Forza L50K offers a massive 8.7 Ncm of torque and a speed range of 1,000-50,000 rpms. The lightweight, ergonomic handpiece weighs only 219 grams ensuring precise control and reduced hand fatigue.





Forza L50K MicroMotor
Order Number: 5021337U0

Motor Speed Range: 1,000-50,000 rpm
 Maximum Torque of 8.7 Ncm
 Brushless
 Weight: 219 grams
 1 Year Warranty



Forza L50K Complete System
Order Number: 5021334U0

System includes Motor, Cord, Control Unit,
 Foot Control and AC Adapter
 Digital Speed Display
 Easy-to-Operate Display Panel
 Foot Pedal or Hand Control Option
 Auto Cruise Function
 1 Year Warranty