



**BRASSELER**  
USA®

BY YOUR SIDE

# ENDOPRO™ Ai1 CORDLESS

The EndoPro Ai1 is an all-in-one endodontic solution. This cordless endodontic handpiece is a balance of all features necessary for a successful RCT. Its ergonomic design with ultra-miniature head will give unprecedented visibility in the posterior.

- ▶ Built in apex locator, eliminating the need for an additional device
- ▶ Pre-programmed to work with all Brasseler files, including reciprocation
- ▶ Compatible with any file on the market
- ▶ Advanced safety features: Apical Action, Auto Start/Stop and a dedicated canal measurement setting



Order No: 5030788U0  
Speed range: 0-3,000 rpm  
Fixed reciprocation setting  
Up to 4.0 Ncm of torque  
Ultra-miniature head size  
1-Year Warranty

*Set Includes: Endo Motor with Contra Angle Head, Charging Base, Measuring Wire, Lip Hook, File Clip, and Probe Lead*

- ▶ Order online at [BrasselerUSA.com](https://www.BrasselerUSA.com) or call **800.841.4522**.  
In Canada call **800.363.3838**.





**BRASSELER**  
USA®

BY YOUR SIDE

# ENDOPRO™ PLUS OBTURATOR GUN

Order No: 5030823U0  
Temp Range: 150-230°C  
Works with BC Pellets™ and other brands  
Weight: 170g  
4-Hour Battery Life  
6 Injection Needle sizes



## EndoPro™ Plus Obturator Gun Injection Needles

Size (Gauge)	Length (mm)	22 mm	24 mm	28 mm
20	22	5030828U0	5030829U0	5030830U0
23	24	—	5030832U0	5030831U0
25	24	—	5030833U0	—

# ENDOPRO™ PLUS DOWN PACK

Order No: 5030797U0  
Cordless Down Pack  
Temp Range: 150-230°C  
Weight: 80g  
4-Hour Battery Life  
9 Tip Sizes



## EndoPro™ Plus Down Pack Tips

Tip Size	Taper	Order No.
35	04	5030799U0
40	04	5030800U0
45	04	5030801U0
50	04	5030802U0
55	06	5030803U0
55	08	5030804U0
55	08L	5030806U0
55	10	5030805U0
60	12	5030807U0

► Order online at [BrasselerUSA.com](http://BrasselerUSA.com) or call **800.841.4522**.  
In Canada call **800.363.3838**.

