

# FORZA V3

PIEZO ULTRASONIC



By Your Side  
in Dentistry




**BRASSER USA**  
SURGICAL POWER & ACCESSORIES

# FORZA V3 PIEZO ULTRASONIC



Dual LED illumination  
delivers clearer vision



Slim design increases  
access and visibility

Elegant and Intuitive. These two words best describe the new Forza V3 piezo ultrasonic from Brasseler USA. At only 33 grams the slim, ergonomic design of the new Forza V3 handpiece provides excellent access while reducing hand fatigue. Stunning LED illumination from dual fiber optics delivers clearer vision and maximum visibility. Equally important is the Piezo Accelerator mode that instinctively provides additional power at the moment you need it most. The compact, space saving unit complements the simple-to-operate interface and can be easily installed in any dental operator. With a wide range of piezo ultrasonic tips available, the Forza V3 is ideal for a variety of clinical procedures including general scaling, perio and endo.

# Elegant. Intuitive. Brasseler.



Piezo accelerator mode  
provides consistent power





Balanced handpiece  
reduces fatigue

## THE PIEZO ADVANTAGE

- *Linear motion provides less sensitivity for the patient.*
- *No vibration in the handpiece increases tactile feel.*
- *Less water required to control the heat.*
- *Wide variety of clinical applications including general scaling, perio, endo, endo surgery and restorative.*
- *Over 80 tips to choose from.*

Actual Size

- Compact Control Unit
- Simple Operation Control Panel
- General, Perio & Endo Applications
- Easy Power & Water Adjustment
- LED Optic and Non-Optic Options



The Control Unit Holder is included with the purchase of every Forza V3 system and provides multiple mounting options. A sliding bracket feature allows the unit to be stored out of way when not being used.



The Piezo Accelerator mode instinctively provides additional power at the moment you need it most. Constantly running in the background, the Piezo Accelerator mode provides consistent power and stable frequency output.

### Forza V3 LED Complete Set

Order No. 5023152U0



### Forza V3 Non-Optic Complete Set


Order No. 5023153U0



- Set Contains:
- Control Unit with Handpiece Cord
  - Handpiece
  - Control Unit Holder
  - Handpiece Holder
  - Water Filter Set
  - Water Supply Connector
  - Foot Control
  - Tip Wrench with Torque Limiter
  - AC Adapter
  - Three Scaling Tips (G1, G9, G11)

## Ultrasonic Tips - Endo

**E4 5011379U0** The short length and non-diamond coated surface makes it ideal for acoustic streaming. Also designed for removing posts in anterior teeth.



**E9D 5012562U0** Diamond coated tip. Multifunctional workhorse excellent for access refinement and troughing.



**E7D 5010237U0** Diamond coated tip. Moderately tapered tip is designed for removing gutta percha and foreign material during retreatment cases and for canal enlargement.




**E14D 5012563U0** Diamond coated tip. The greater taper and shorter length is designed for troughing and access refinement in smaller canals.



**E9 5012561U0** Effective for removing posts in molars and premolars and for acoustic streaming.




**E15D 5019344U0** The diamond coated head design is slightly smaller than a #2 round bur (.095). Effective for troughing and pulp stone removal.




## Ultrasonic Tips - Perio/Root Planing


**PI 5012062U0** Effective for removal of light calculus in tight difficult to reach interproximal areas.




**P20 5010242U0** Effective for removal of light deposits and biofilm. Extended, thin tip designed for deep pockets.




**P1D 5012063U0** Diamond coated tip. Effective for removal of light calculus in tight difficult to reach interproximal areas.




**P26R 5022654U0** Effective for removal of light to medium deposits and biofilm. Bladed tip and curved shank design adapts to the root anatomy.



**P10 5010241U0** Effective for difficult to remove flinty, burnished calculus. Bladed tip resembles an anterior Universal sickle.



**P26L 5022655U0** Effective for removal of light to medium deposits and biofilm. Bladed tip and curved shank design adapts to the root anatomy.



# Ultrasonic Tips - Scaling

**G1 5010218U0** Effective for removal of medium to heavy deposits.



**G4 5010221U0** Effective for removal of light to medium deposits and biofilm.




**G6 5010223U0** Effective for removal of light deposits and biofilm during periodontal maintenance.




**G8 5020173U0** Effective for heavy deposit removal.



**G9 5022653U0** Effective for removal of light to medium deposits and biofilm. Long thin straight tip designed for access in deep pockets.



**G11 5010225U0** Effective for removal of medium deposits.



## Forza V3 LED Handpiece Only

Order No. 5024337U0



## Forza V3 Non-Optic Handpiece Only

Order No. 5024338U0



## E-Tip Wrench



Order No. 5010484U0

## CR-10 Torque Limiting Wrench



Order No. 5019770U0

By Your Side  
in Dentistry



**BRASSLER USA**  
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.  
In Canada call 800.363.3838 or fax 866.330.4454.  
Visit our website at [BrasselerUSA.com](http://BrasselerUSA.com)

