

ANOTHER INNOVATION FROM BRASSELER USA...

Dia^{comp}TM *FEATHER LITE*TM Composite

Two-Step Composite Polishing Intra-Oral System

NEW



FLDC5M
Step 1:
Medium Green



FLDC5
Step 2:
Fine Gray

TIP: Following the final contour with ET[®] carbides or diamonds, use the Dia^{comp}TM Feather LiteTM two-step polishers at slower speeds and with water to polish composites.

Ronald E. Goldstein, DDS

By Your Side
in Dentistry



BRASSELER USA[®]
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.
In Canada call 800.363.3838 or fax 866.330.4454.
Visit our website at BrasselerUSA.com



Dia compTM FEATHER LITETM Composite

Two-Step Composite Polishing Intra-Oral System

TWO STEPS LEAVE A
HIGH SHINE ON ALL
COMPOSITES

EASILY ADAPTS
TO ANY
SURFACE

IDEAL FOR
OCCLUSAL
SURFACES

CREATES AN
ENAMEL-LIKE FINISH
ON COMPOSITES

PRESERVES SURFACE
DIFFERENTIATION
WHILE POLISHING

FLEXIBLE
POLISHING
SPIRALS

STERILIZABLE
AND
RE-USABLE

DIAMOND
IMPREGNATED



Order No.	FLDC5M	FLDC5	FLDC5-5M
Head Size (1/10 mm)	140	140	140
Cutting Length (mm)	1.6	1.6	1.6
Grit	medium	fine	medium / fine

21 - RA - Right-Angle Latch (22 mm)

Qty/Pack	3	3	6 (3 each)
----------	---	---	------------

Sterilization: Autoclave

Rec: 5-8,000 RPM*

*for final polish reduce speed by half (2-4,000 RPM)

STEP 1



STEP 2



Photos courtesy of
Ronald E. Goldstein, DDS

By Your Side
in Dentistry



BRASSELER USA
DENTAL INSTRUMENTATION



To order call 800.841.4522 or fax 888.610.1937.
In Canada call 800.363.3838 or fax 866.330.4454.
Visit our website at BrasselerUSA.com

