



**BRASSELER**  
USA®


**Versa***flex*™

## ► A New Spin on Snap-On Polishers

Brasseler USA introduces an innovative new spin on basic, plastic snap-on polishers with the new VersaFlex™ diamond-coated disposable polishing system. Unlike traditional mylar (plastic) polishers, this versatile polishing system is outfitted with a durable diamond coating that delivers unprecedented results. **It easily and effectively finishes and polishes a wide variety of restorative materials, including composites, all-ceramics and metal.**



The flexible design of the VersaFlex discs easily conforms to the individual contours of each tooth, providing unrestricted access to even hard-to-reach locations. The addition of the corresponding mylar strips will help you gain access to tight interproximal areas. Available direct only from Brasseler USA, let VersaFlex enhance your restorations today!

- Four Grit Options: Super-Coarse, Coarse, Medium and Fine
- Two Discs Sizes: 12mm and 10mm
- Two 2-Grit Mylar Strips: Super-Coarse/Coarse and Medium/Fine
- Two Snap-On Mandrel Options: Right-Angle Latch (ra) and Straight Handpiece (hp)



Photographs courtesy of Douglas A. Terry, DDS

**BY YOUR SIDE**






## K0315

- 25 VFD1LG Large 12mm Super-Coarse Green Discs
- 25 VFD1SM Small 10mm Super-Coarse Green Discs
- 25 VFD2LG Large 12mm Coarse Blue Discs
- 25 VFD2SM Small 10mm Coarse Blue Discs
- 25 VFD3LG Large 12mm Medium Pink Discs
- 25 VFD3SM Small 10mm Medium Pink Discs
- 25 VFD4LG Large 12mm Fine Gray Discs
- 25 VFD4SM Small 10mm Fine Gray Discs
- 25 VFS1W Super-Coarse Green/Coarse Blue Strips
- 25 VFS2W Medium Pink/Fine Gray Strips
- 1 310.21 Right-Angle Latch Snap-On Mandrel
- 1 310.11 Straight Handpiece Snap-On Mandrel



Head Size ( $\frac{1}{10}$ mm)	12.0	10.0	12.0	10.0	12.0	10.0	12.0	10.0
Thickness (mm)	0.21	0.21	0.17	0.17	0.14	0.14	0.11	0.11
Grit	Super-Coarse	Super-Coarse	Coarse	Coarse	Medium	Medium	Fine	Fine
Qty/Pk	100	100	100	100	100	100	100	100
Unmounted	VFD1LG	VFD1SM	VFD2LG	VFD2SM	VFD3LG	VFD3SM	VFD4LG	VFD4SM



	Width (mm)	Length (mm)	Grit		Thickness (mm)		Qty/Pk
VFS1W	4	150	Super-Coarse	Coarse	0.21	0.17	100
VFS2W	4	150	Medium	Fine	0.14	0.11	100



## 310

Length (mm)	22.0	44.5
Qty/Pk	5	5
21-RA-Right Angle Latch	310.21	—
11- HP- Straight Handpiece	—	310.11

Discs and strips are single patient.

Recommended disc speed 7,000 - 10,000 rpms.

Visit our website at [BrasselerUSA.com](http://BrasselerUSA.com)

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838.

