



BRASELER
USA®

NEW INSTRUMENTATION ORGANIZATION & STORAGE

BRASELER USA HAND INSTRUMENT CASSETTES




PROCEDURAL CASSETTES

Procedural cassettes help to protect, control, and organize inventory using a systematic approach to help your staff achieve maximum efficiency in patient care.


Storing instruments by procedure type maximizes efficiency, reduces reprocessing time, and protects instruments from breakage. Also, reprocessing instruments inside a closed cassette minimizes staff exposure to sharp instrument tips.

Manufactured using highest grade, passivated stainless steel to maintain aseptic appearance after repeated sterilization cycles



*Not available on all cassettes, see cassette descriptions for features

SEVEN SIZES TO MEET YOUR OFFICE'S NEEDS




5 Count Cassette

Order number: HINSTC05R

Perfect size for exam and routine hygiene maintenance setups

- Single hinge lid
- Lid Clip included
- Accessory compartment
- Color coding tabs
- Dimensions 8" X 3 5/8" X 1 1/4"




Handpiece Cassette with Mat

Order number: HINSTC05

Holds two handpieces, or can be used as a small prophy or exam cassette holding up to 5 instruments

- Single hinge lid
- No interior instrument rack
- Holds 5 or more hand instruments, depending on size
- Dimensions 8 X 2 1/2" X 11/16"




7 Count Slimline Cassette

Order number: HINSTC07S

Reduced profile ideal for small sterilizers

- Single hinge lid
- No interior instrument rack,
- Holds 7 or more hand instruments, depending on handle size
- Dimensions 8 X 2 1/2" X 11/16"




7 Count Cassette

Order number: HINSTC07

The ideal rack size for smaller procedure setups

- Double hinge lid
- Lid Clip included
- Accessory compartment
- Color coding tabs
- Dimensions 8" X 4 3/4" X 1 1/4"




9 Count Cassette

Order number: HINSTC09

Designed for expanded procedure setups, like SRP, composite, etc.

- Double hinge lid
- Needle re-capper
- Accessory compartment
- Color coding tabs
- Dimensions 8 X 6 X 1 ¼"




10 Count Cassette

Order number: HINSTC10

Holds 10 hand instruments and features side compartment for storage of syringes, sharpening stones, and ring-handled or odd-sized instruments.

- Double hinge lid
- Needle re-capper
- Side compartment size 3 1/8" X 7 15/16"
- Color coding tabs
- Dimensions 8 X 9 X 1 ¼"



16 Count Cassette

Order number: : HINSTC16

Generous size is just right for specialty setups and features side compartment for syringes, sharpening stones, and ring-handled or odd-sized instruments.

- Double hinge lid
- Needle re-capper
- Size compartment size 3 ½" X 8"
- Color coding tabs
- Dimensions 8 3/16" X 10 ¾" X 1 3/16"

Code Rings

Color code instrument rings available in small/regular (1/8 inch) and large (1/4 inch). Single color packages contain 50 pieces. Assorted color packages contain five (5) pieces of each color, a total of 50 pieces.



Code Rings Color	Product Order Number Small Rings	Product Order Number Large Rings
Blue	CR-BLS	CR-BLL
Pink	CR-PKS	CR-PKL
White	CR-WS	CR-WL
Black	CR-BKS	CR-BKL
Red	CR-RS	CR-RL
Yellow	CR-YS	CR-YL
Green	CR-GRNS	CR-GRNL
Orange	CR-OS	CR-OL
Purple	CR-PUS	CR-PUL
Gray	CR-GRYS	CR-GRYL
Assortment	CR-ASSTS	CR-ASSTL



BRASSLER
USA®



B-4742-US-XXX-11.16-XXX