



BRIOSHINE FEATHER LITE

SINGLE-USE UNIVERSAL COMPOSITE POLISHING SYSTEM

- Adapts to most tooth surfaces
- Reduces time spent polishing
- Fewer instruments in tray set-ups
- Recommended for use at slow speed (4-6,000 rpm) and light water



UNIVERSAL COMPOSITE STEP 1

Pre-polishing composite restorations



UNIVERSAL POLISHING STEP 2

Final, universal polisher that provides an excellent finish to all restorations



	CON	E.S	
Order No.	FLBC0-1	FLBF-2	
Head Size ($^{1}/_{10}$ mm)	100	100	
21 - RA - Right-Angle Latch (22 mm)			
Qty/Pack	30	30	

Order online at BrasselerUSA.com or call 800.841.4522. In Canada call 800.363.3838.



©2018 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-4967-01.18-20M-PRE





BRIOSHINE FEATHER LITE

SINGLE-USE UNIVERSAL CERAMIC POLISHING SYSTEM

- Adapts to most tooth surfaces
- Reduces time spent polishing
- Fewer instruments in tray set-ups
- Recommended for use at slow speed (4-6,000 rpm) and light water



UNIVERSAL CERAMIC STEP 1

Pre-polishing zirconia, lithium disilicate and porcelain/PFM restorations



UNIVERSAL POLISHING STEP 2

Final, universal polisher that provides an excellent finish to all restorations



Head Size (1/10 mm)	100	100		
21 - RA - Right-Angle Latch (22 mm)				
Qty/Pack	30	30		

Order online at BrasselerUSA.com or call 800.841.4522. In Canada call 800.363.3838.



©2018 Brasseler USA. All rights reserved. Printed in the U.S.A.