



**BRASSELER**  
USA®

# BRASSELER USA ENDO SOLUTIONS



BY YOUR SIDE

# Endo SOLUTIONS

Studies have shown that 15% of store bought NaOCl contains less than 1% available chlorine (which is needed for tissue dissolution and disinfection).\* Brasseler USA's Endo Solutions adhere to strict FDA & ISO quality control measures.

## 5% SODIUM HYPOCHLORITE (16 oz)

Retail: \$39.95

Item Number: 5026552U0

- Eliminates bacteria
- Used to irrigate during and after root canal debridement

## 2% CHLORHEXIDINE GLUCONATE (16 oz)

Retail: \$52.95

Item Number: 5026553U0

- Aids in post-op sensitivity
- Increases bond strength

## 17% EDTA (16 oz)

Retail: \$52.95

Item Number: 5026554U0

- Opens dentinal tubules
- Prepares dentin walls for obturation
- Removes smear layer



Exclusive luer-lock cap reduces cross-contamination and improves solution efficacy.



## IRRIGATION SYRINGES

Brasseler USA's irrigation syringes are ideal for surgical site irrigation and compatible with all luer-lock style needles and hubs. Disposable, non-sterile. Not for injection.

Retail: \$14.75

3 CC SYRINGE (50/Pk)

Item Number: 5026555U0

6 CC SYRINGE (50/PK)

Item Number: 5026882U0

## IRRIGATION NEEDLES

Designed with a closed-end cannula and a side window dispersal point, Brasseler USA's Irrigation Needles offer a unique design for unmatched effectiveness, patient comfort and safety in irrigating root canals.

Retail: \$22.25

27 GAUGE IRRIGATION NEEDLES (20/PK)

21mm (yellow)

Item Number: 5026557U0

27mm (green)

Item Number: 5026558U0


31 GAUGE IRRIGATION NEEDLES (20/Pk)

21mm (yellow)

Item Number: 5026680U0

27mm (green)

Item Number: 5026681U0



\*Quintessence Intl 2014;45:467-474; doi:10.3290/j.qi.a31805)

# ENDO SOLUTIONS CADDY

Retail: \$102.95

Item Number: 5026556U0

Dimensions: 12.5"L x 10"W x 7.25"H

Brasseler USA's Endo Solutions Caddy is the smart way to store your irrigation and endo supplies. The caddy features three forward-angled slots for bottles, allowing you to take full advantage of the luer-lock cap. The system virtually eliminates waste and the risk of exposing solutions to damaging air, light and heat. The caddy also features three large base compartments which can store various endo supplies including tips, files, gutta percha and BC Sealer. The top of the caddy contains slots to accommodate up to twelve syringes (8 x 3cc & 4 x 6cc).



## PULP VITALITY SPRAY (6 oz)

Retail: \$23.50

Item Number: 5026551U0

- A quick and effective way to check cold sensitivity
- Significantly colder than ethyl chloride or ice
- For accuracy, use with micro applicators or cotton pellets
- Directional nozzle for easy placement on micro applicator or cotton pellet







DIAMONDS



CARBIDES



POLISHERS



PROCEDURE SYSTEMS



BIOCERAMICS



Since 1976, our absolute focus has been to develop products and provide services that support the practice of core operator procedures. Today, Brasseler USA® is the premier dental instrumentation company in North America offering the most extensive selection of dental and surgical instrumentation under one brand.



HANDPIECES



ENDODONTICS



SMALL EQUIPMENT



HYGIENE



HAND INSTRUMENTS



ORAL SURGERY



LABORATORY INSTRUMENTS



Visit our website at [BrasselerUSA.com](http://BrasselerUSA.com)

To order call 800.841.4522 or fax 888.610.1937. In Canada call 800.363.3838

