



**BRASSELER  
USA®**

BY YOUR SIDE

# True™ Dental Dam

Isolate With Confidence.

True Dental Dam provides the fit, seal and isolation you demand for better clinical outcomes.

## ► Simple Placement

- Easily punched, stretched, and placed, supporting effective and efficient isolation set-up procedures.
- Smoothly glides through contacts and tucks into the sulcus, allowing for true placement confidence.

## ► Durability & Protection

- Exceptionally durable and tear resistant, delivering the seal and tissue protection you expect and help to reduce contamination to create a clean, clear operative field.
- Available in blue (latex) and purple (latex-free) to provide better color contrast with dentition.

► Order online at [BrasselerUSA.com](https://www.BrasselerUSA.com) or call **800.841.4522**.



©2021 Brasseler USA. All rights reserved. Printed in the U.S.A.



B-5537-08.21

# True™ Dental Dam

## ▶ Product Offering

With both latex and latex-free options available in two sizes and two gauges, there is a True Dental Dam sure to suit your needs.



Order No.	Size (in.)	Gauge	Material Type	Color	Sheets Per Pack	Thickness (mm)
838252-B	5x5	Medium	Latex	Blue	52	0.23
838236-B	6x6	Medium	Latex	Blue	36	0.23
838200-B	5x5	Medium	Latex-Free	Purple	50	0.22
838201-B	6x6	Medium	Latex-Free	Purple	50	0.22
838253-B	5x5	Heavy	Latex	Blue	52	0.28
838237-B	6x6	Heavy	Latex	Blue	36	0.28
838202-B	5x5	Heavy	Latex-Free	Purple	50	0.25
838203-B	6x6	Heavy	Latex-Free	Purple	50	0.25

## True™ Dental Dam & EndoSequence® Flowable Dam Material

Excellent Isolation, Together.

Use True Dental Dam for primary isolation and EndoSequence Flowable Dam Material for secondary isolation to help reduce gaps and improve the seal.

### ▶ EndoSequence Flowable Dam Material

- Flowable yet stackable consistency
- Light cures quickly and stays cool to the touch
- Flexible and tear-resistant
- Stays put while remaining easy to remove
- Can be applied to gingiva or tooth surfaces before root canal therapy, bleaching, sandblasting, or other procedures requiring intra-oral protection



EndoSequence Flowable Dam Material  
4 3mL Syringes  
40 Tips

Order No. 5028289U0

*Clinical image courtesy of Allen Ali Nasseh, DDS, MMSc.*

▶ Order online at [BrasselerUSA.com](http://BrasselerUSA.com) or call **800.841.4522**.

