
BRASSELER USA®

BY YOUR SIDE

SOFT TISSUE TRIMMERS

The new Soft Tissue Trimmer from Brasseler USA features a ceramic tip for trimming gingival tissue with minimal bleeding. The Soft Tissue Trimmer can be used in a wide range of soft tissue procedures and is an effective alternative to scalpels or electro-surgery.



Order No.	STT	STTL
Head Size (1/10 mm)	016	016
Cutting Length (mm)	5.0	8.0
31 - FG - Friction Grip (19 mm)		
Qty/Pack	1	1

Recommended speed:
300,000 – 500,000 rpm
in high speed air handpiece

Use without water spray


Light hand pressure



Range of Indications:

- Widening the sulcus for impressions
- Exposing deep cervical cavities
- Exposing impacted teeth
- Exposing intra-osseous implants
- Crown Lengthening
- Papillectomies
- Gingiva contouring

Photo courtesy of Jeff Lineberry, DDS, FAGD, FECOI, AAACD



▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2019 Brasseler USA. All rights reserved. Printed in the U.S.A.