

solutions™


# Dr. Ronald E. Goldstein ET® Carbide Combo Kit




## 4200



| US No. | ET3 | ET3F | ET3UF | ET4 | ET4F | ET4UF | ET6 | ET6F | ET6UF | ET9 |
|---|---|---|---|---|---|---|---|---|---|---|
| Figure No. | H132 | H132F | H132UF | H133 | H133F | H133UF | H134 | H134F | H134UF | H135 |
| Head Size (1/10 mm) | 008 | 008 | 008 | 010 | 010 | 010 | 014 | 014 | 014 | 014 |
| Cutting Length (mm) | 3.0 | 3.0 | 3.0 | 4.0 | 4.0 | 4.0 | 6.0 | 6.0 | 6.0 | 9.0 |



| US No. | ET9 | ET9UF | 7408/0S1 | 0S1F | 0S1UF | 7801/0S2 | 0S2F | 0S2UF | 0S3 | 0S3F | 0S3UF |
|---|---|---|---|---|---|---|---|---|---|---|---|
| Figure No. | H135F | H135UF | H379 | H379F | H379UF | H247 | H247F | H247UF | H247 | H247F | H247UF |
| Head Size (1/10 mm) | 014 | 014 | 023 | 023 | 023 | 009 | 009 | 009 | 007 | 007 | 007 |
| Cutting Length (mm) | 9.0 | 9.0 | 4.2 | 4.2 | 4.2 | 3.2 | 3.2 | 3.2 | 3.2 | 3.2 | 3.2 |


Order online at BrasselerUSA.com
US: 800.841.4522 or fax 888.610.1937
Canada: 800.363.3838 or fax 866.330.4454

BRASSELER USA®

©2021 Brasseler USA. All rights reserved.

B-3732B-REV2-5022192L0-05.21

# solutions™ 4200 ET® Dr. Ronald E. Goldstein Carbide Combo Kit

| | Order No. | Head Size | Grit | Usage |
|---|---|---|---|---|
|  | ET3 | 008 | F | Large amount gingival/interproximal finishing of direct bonded restorations |
|  | ET3F | 008 | EF | Minimal amount of gingival or proximal finishing of direct resin restorations |
|  | ET3UF | 008 | UF | Final subgingival/interproximal finishing of bonded/provisional/ceramic restorations |
|  | ET4 | 010 | F | Gross gingival finishing of direct bonded restorations |
|  | ET4F | 010 | EF | Minimal gingival finishing of direct bonded or provisional restorations is required |
|  | ET4UF | 010 | UF | Final finishing of direct bonded, provisional, or ceramic restorations |
|  | ET6 | 014 | F | Gross labial/lingual finishing of medium sized teeth |
|  | ET6F | 014 | EF | Labial/lingual minimal finishing anterior/posterior teeth |
|  | ET6UF | 014 | UF | Final finishing or placing texture on labial surfaces of bonded or provisional restorations |
|  | ET9 | 014 | F | Gross contouring of labial surfaces of direct resin or provisional restorations |
|  | ET9F | 014 | EF | Minimal contouring or finishing of labial or lingual surfaces on longer teeth |
|  | ET9UF | 014 | UF | Texturing/final finish of labial surfaces on longer bonded/provisional restorations |

| | Order No. | Head Size | Grit | Usage |
|---|---|---|---|---|
|  | 7408/OS1 | 023 | F | Gross lingual/occlusal finishing and occlusal equilibration procedures |
|  | OS1F | 023 | EF | Minimal finishing lingual/occlusal surfaces of bonded/provisional restorations or slight adjustment of natural teeth |
|  | OS1UF | 023 | UF | Final occlusal/lingual finishing bonded restorations/occlusal adjustments temporary provisional materials |
|  | 7801/OS2 | 009 | F | Occlusal anatomy in posterior bonded or provisional restorations |
|  | OS2F | 009 | EF | Initial placement of occlusal anatomy when lesser amounts are required |
|  | OS2UF | 009 | UF | Final polishing/finishing in posterior bonded/temporary restorations |
|  | OS3 | 007 | F | Supplemental occlusal anatomy in posterior bonded restorations/provisional restorations |
|  | OS3F | 007 | EF | Placement of supplemental occlusal anatomy for fine grooves |
|  | OS3UF | 007 | UF | Finish supplemental occlusal anatomy in bonded/temporary restorations |

BRASSELER USA®

Brasseler U.S.A. Dental, LLC
One Brasseler Boulevard
Savannah, GA 31419 U.S.A.