BRASSELER USA®

BY YOUR SIDE

IMPLANTPRO®
TITANIUM PROBES

Ideal for use on both implants and natural dentition

Designed with the ultimate combination of innovative design and superior ergonomics, ImplantPro® Titanium Probes feature precise color markings on a completely smooth surface for a more gentle diagnostic evaluation and improved patient comfort.

- Bright zirconium nitride coating provides superior color contrast against black millimeter markings for easy reading.
- All feature Brasseler USA's distinctive, lightweight PEEK (Polyether Ether Ketone) handles for outstanding ergonomics.

TICPOW
Michigan "O" with Williams markings



TICPUNC15
UNC15 Probe



TICP12
Marquis Probe



Contact your Brasseler USA representative today to learn more about the ImplantPro® product line.

▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2017 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-4816-04.16-30M-RRD

BRASSELER USA®

BY YOUR SIDE

IMPLANTPRO®
TITANIUM 204S SCALER AND LTS CURETTES

Providing the ideal balance of safety and efficiency

ImplantPro® Titanium Scalers and Curettes are made from non heat-treated 6Al-4V titanium, keeping the Rockwell C hardness at a low 25-31 HRC, comparable to most implant abutments. All feature Brasseler USA's distinctive, lightweight PEEK (Polyether Ether Ketone) handles for outstanding ergonomics. Special blade finishing with long-lasting, polished edges ensures safe instrumentation around implants. Due to this unique blade finishing, hand sharpening is not recommended.

204S SCALER

· Features thinner blades with less curvature than a traditional 204S

TI204S
Titanium 204S



LTS UNIVERSAL CURETTES

· Features elongated terminal shanks and thinner, shorter blades to reduce tissue distension and ease subgingival insertion.

TI4R/4L
Titanium 4R/4L



TI1/2
Titanium 1/2



TI11/12
Titanium 11/12



TI13/14
Titanium 13/14



Contact your Brasseler USA representative today to learn more about the ImplantPro® product line.

▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2017 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-4816-04.16-30M-RRD