BRASSELER USA®

BY YOUR SIDE


NEW!

C SERIES™

TRIMMING & FINISHING CARBIDES

- Durable, single-piece* design
- New blade geometry
- Ideal for shaping, contouring and finishing of all composites

**Limited to 016 head size and smaller*



Order No.	7404CM	7406CM	7408CM	7901CM	7902CM	7903CM
Head Size (1/10 mm)	014	018	023	009	010	012
Cutting Length (mm)	3.2	3.5	3.6	3.8	4.1	4.0
Friction Grip						
31 - FG (19mm)	FG	FG	FG	FG	FG	FG



- Order online at **BrasselerUSA.com** or call **800.841.4522**. In Canada call **800.363.3838**.

©2019 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-5273-08.19