
BRASSELER USA®

BY YOUR SIDE

ZIRPRO™ ALL CERAMIC FINISHING INSTRUMENTS



- Diamond-impregnated and ceramic-bonded
- For adjusting, contouring, and finishing of modern ceramic materials
- Recommended for use on Zirconia in the sintered state and on IPS e.max® in the pre-crystallized (blue) state
- Variety of shapes to adjust and finish all surfaces of the restoration
- Cool and efficient finishing



Quick reduction of sprue at 10,000 rpm

Adjusting and contouring at 5,000 rpm



Order No.	4000.HP	4002.HP	4008.HP	4010.HP	4011.HP	4012.HP
Head Size (1/10 mm)	220	220	035	115	140	065
Cutting Length (mm)	4.5	2.0	11.5	3.0	9.0	3.5
Grit	Coarse	Coarse	Coarse	Coarse	Coarse	Coarse
11 - HP - Straight Handpiece (44.5 mm)						
Qty/Pack	1	1	1	1	1	1

Recommended Speed: 10,000 – 15,000 rpm; Max Speed: 35,000 rpm

Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2019 Brasseler USA. All rights reserved. Printed in the U.S.A.
IPS e.max® is a registered trademark of Ivoclar Vivadent®