BRASSELER USA®

BY YOUR SIDE

# BRIO·LS

*Same Warranty, Same Performance, Significant Savings.*

Brasseler USA now offers an affordable alternative to expensive replacement heads with a line of compatible Low Speed Heads for your Midwest®, KaVo, and Star® Low Speed Attachments. The Brio LS Series seamlessly integrates into your existing Low Speed Air set ups.

## MIDWEST® LOW SPEED ATTACHMENT COMPATIBLE




**Brio LS-MW-L**
**Order Number: 5028352U0**

For Latch Type Burs
Push Button Head
1-Year Warranty




**Brio LS-MW-F**
**Order Number: 5028353U0**


For Friction Grip Burs
Push Button Head
1-Year Warranty

## KAVO® LOW SPEED ATTACHMENT COMPATIBLE




**Brio LS-KV-L**
**Order Number: 5028348U0**

For Latch Type Burs
Push Button Head
1-Year Warranty




**Brio LS-KV-F**
**Order Number: 5028349U0**

For Friction Grip Burs
Push Button Head
1-Year Warranty

## STAR® MOTOR-TO-ANGLE ADAPTER COMPATIBLE




**Brio LS-ST-L**
**Order Number: 5028350U0**

For Latch Type Burs
Push Button Head
1-Year Warranty




**Brio LS-ST-F**
**Order Number: 5028351U0**

For Friction Grip Burs
Push Button Head
1-Year Warranty

*See reverse side for compatibility chart.

*Midwest® is a registered trademark of DENTSPLY International. KaVo® is a registered trademark of Kaltenbach & Voigt GmbH, Germany. Star® is a registered trademark of DentalEZ, Inc.*

Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2020 Brasseler USA. All rights reserved. Printed in the U.S.A.

# BRIO™ LS

## BRIO LS COMPATIBILITY CHART

| Item Number | Model |
|---|---|
| 5028352U0 | BRIO LS-MW-L LATCH HEAD (for latch type burs) |
| Replaces | Midwest® Push-Button Latch Type |
| For Use With | Midwest® Contra Angle Sheath<br>Henry Schein® Maxima®-M Contra Angle Sheath |

| Item Number | Model |
|---|---|
| 5028353U0 | BRIO LS-MW-F FG HEAD (for friction grip burs) |
| Replaces | Midwest® Push-Button Friction Grip |
| For Use With | Midwest® Contra Angle Sheath<br>Henry Schein® Maxima®-M Contra Angle Sheath |

| Item Number | Model |
|---|---|
| 5028348U0 | BRIO LS-KV-L LATCH HEAD (for latch type burs) |
| Replaces | KaVo® L68 B & 68LU Intra Lux™ Attachment Head |
| For Use With | KaVo® Mastermatic™ attachments M20L, M07L, M29L.<br>KaVo® Intraflex Lux™ 2 2307LN, 2320LN, 2325LN motors<br>KaVo® Gentlepower™ 20LP, 10LP & 7LP attachments<br>Henry Schein® Maxima®-K Contra Angle Sheath |

| Item Number | Model |
|---|---|
| 5028349U0 | BRIO LS-KV-F FG HEAD (for friction grip burs) |
| Replaces | KaVo® L80 & 80LD Intra Lux™ Attachment Head |
| For Use With | KaVo® Mastermatic™ attachments M20L, M07L, M29L.<br>KaVo® Intraflex Lux™ 2 2307LN, 2320LN, 2325LN motors<br>KaVo® Gentlepower™ 20LP, 10LP & 7LP attachments<br>Henry Schein® Maxima®-K Contra Angle Sheath |

| Item Number | Model |
|---|---|
| 5028350U0 | BRIO LS-ST-L LATCH HEAD WITH ELBOW (for latch type burs) |
| Replaces | Star® Lube Free Auto Latch Angle<br>Star® Ball Bearing Auto Latch Angle |
| For Use With | Star® Motor-To-Angle Adapter<br>Henry Schein® Maxima®-S Contra Angle Adapter |

| Item Number | Model |
|---|---|
| 5028351U0 | BRIO LS-ST-F FG HEAD WITH ELBOW (for friction grip burs) |
| Replaces | Star® Lube Free Auto Latch Angle<br>Star® Ball Bearing Auto Latch Angle |
| For Use With | Star® Motor-To-Angle Adapter<br>Henry Schein® Maxima®-S Contra Angle Adapter |

HANDPIECE REPAIR X

## DID YOU KNOW?

- Brasseler repairs all major models of high speed, low speed, and electric attachments.
- Brasseler now offers replacement low speed heads for Midwest®, KaVo®, and Star® low speeds.
- Replacement low speed heads come with a warranty EQUAL TO OR GREATER THAN the original manufacturer's warranty!



*Not sure if we can repair your handpiece?*
*Contact your sales representative or contact Brasseler USA directly!*
*1.800.841.4522 (ext. 4675)*

*Midwest® is a registered trademark of DENTSPLY International. KaVo®, Intra Lux™, Mastermatic™, and Gentlepower™ are trademarks of Kaltenbach & Voigt GmbH, Germany. Star® is a registered trademark of DentalEZ, Inc. Henry Schein® and Maxima® are registered trademarks of Henry Schein, Inc.*

BRASSELER USA®

©2020 Brasseler USA. All rights reserved. Printed in the U.S.A.