
BRASSELER USA®

BY YOUR SIDE

DURACUT™ DIAMONDS

Using less efficient diamonds can make tooth preparation, endodontic access, and failed zirconia crown removal time-consuming, increasing chair time. Brasseler's specialized DuraCut Diamonds have improved efficiency and effectively cut modern ceramics and natural dentition with a smoother and longer-lasting cutting performance.






▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2021 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-5565-12.21 REV 0

DuraCut™ Diamonds



31 - FG

Product Offering

DuraCut Diamonds are offered in a large selection of the most popular shapes and sizes for multiple procedures. All DuraCut Diamonds are available in convenient 5-packs.

ENDO ACCESS




Head Size (1/10 mm)	012	014	018	012	014	018	023
Carbide Equivalent	–	–	–	#3	#4	#6	#8
Cutting Length (mm)	10.0	10.0	10.0	–	–	–	–
Friction Grip							
31 - FG (19 mm)							
DuraCut™ (DC)	15802DC	15802DC	15802DC	6801DC	6801DC	6801DC	6801DC

BREAKING CONTACT / OCCLUSAL/LINGUAL REDUCTION / SHOULDER / MODIFIED SHOULDER



Head Size (1/10 mm)	014	014	040	033	023	023	014	016	016	018
Cutting Length (mm)	8.0	10.0	1.8	4.5	4.2	5.0	8.0	8.0	8.0	8.0
Friction Grip										
31 - FG (19 mm)										
DuraCut™ (DC)	6858DC	6850DC	6909DC	6811DC	6379DC	6368DC	6847DC	6847DC	6847KRDC	6847KRDC

CHAMFER / MODIFIED CHAMFER / FLAME



Head Size (1/10 mm)	014	016	018	018	016	014	016	018	014
Cutting Length (mm)	8.0	8.0	8.0	9.0	10.0	8.0	8.0	8.0	8.0
Friction Grip									
31 - FG (19 mm)									
DuraCut™ (DC)	6856DC	6856DC	6856DC	6856LDC	6850DC	6878KDC	6878KDC	6878KDC	6862DC

Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

BRASSELER USA®

©2021 Brasseler USA. All rights reserved. Distributed by Brasseler U.S.A. Dental, LLC.Printed in the U.S.A.