

BRASSELER USA®

BY YOUR SIDE

DIALITE® DIAMONDS

Multi-layered Dialite Diamonds are designed for adjusting and removing ceramic restorations, including zirconia, lithium disilicate, and feldspathic porcelain. Adjust restorations atraumatically with Dialite Diamonds to minimize flaw induction and potential ceramic fracture. Remove failed ceramic restorations efficiently and comfortably. Use them in the dental laboratory to shape and create surface texture in pre-sintered (green state) zirconia and pre-crystalized (blue block) lithium disilicate.






▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2022 Brasseler USA. All rights reserved. Printed in the U.S.A.


Dialite® Diamonds



31 - FG

Product Offering

Dialite Diamonds can be used for ceramic restoration adjustment and removal and are available in convenient 5-packs.
Recommended speeds: adjusting—40,000-60,000 rpm; removal—100,000 rpm



369D

Head Size (1/10 mm)	025
Cutting Length (mm)	5.5
Friction Grip	
31 - FG (19 mm)	
Fine	8369DF
Extra-Fine	369DEF



379D

Head Size (1/10 mm)	023
Cutting Length (mm)	4.2
Friction Grip	
31 - FG (19 mm)	
Fine	8379DF



801LD

Head Size (1/10 mm)	014
Carbide Equivalent	#4
Friction Grip	
31 - FG (19 mm)	
Fine	8801LDF
Extra-Fine	801LDEF



856D

Head Size (1/10 mm)	016
Cutting Length (mm)	8.0
Friction Grip	
31 - FG (19 mm)	
Fine	8856DF
Extra-Fine	856DEF




KS1D

Head Size (1/10 mm)	012
Cutting Length (mm)	7.0
Friction Grip	
31 - FG (19 mm)	
Fine	KS1DF



274D

Head Size (1/10 mm)	016
Cutting Length (mm)	3.5
Friction Grip	
31 - FG (19 mm)	
Fine	8274DF



USE WITH DIALITE® DIAMONDS

Converter Mandrel (10/Pack)
311FGHP

The FG to HP Converter Mandrel allows laboratory technicians to use friction grip (FG) instruments in devices designed for straight handpiece (HP) instruments, providing access to a wider variety of instrumentation without two different handpiece units.

Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2022 Brasseler USA. All rights reserved. Distributed by Brasseler U.S.A. Dental, LLC.Printed in the U.S.A.