BRASSELER USA®

BY YOUR SIDE

BRIO 550

REDUCE AEROSOLS, NOT PERFORMANCE

Specifically engineered for effective aerosol reduction, the powerful Brio550 high speed air handpiece line helps reduce airborne particle emissions by an average of 90%* compared to typical high speed handpieces.



- INCREASED POWER, UP TO 21 WATTS
- WEAR-RESISTANT COATING
- AEROSOL REDUCTION TECHNOLOGY
- ONE-YEAR WARRANTY


*Average reduction based on tests conducted in an aerosol testing chamber and measured by laser-light scattering photometry.

- Order online at **BrasselerUSA.com** or call **800.841.4522**.
  In Canada call **800.363.3838**.

BRIO 550


AEROSOL REDUCTION

Reduces emission of airborne micro-droplets by an average of more than 90% compared to typical highspeed air handpieces.



ANTI-RETRACTION

Active sealing system prevents suck-back of fluid and particles to reduce risk of cross-contamination of air and water lines.





Brio550L™ Fiber Optic

Order Number: 5029832U0

For Brasseler/NSK® Connections
21 Watts of Power
1-Year Warranty

Brio550™ Non-Optic

Order Number: 5029833U0

For Brasseler/NSK® Connections
21 Watts of Power
1-Year Warranty

Brio550KL™ Fiber Optic

Order Number: 5029834U0

For KaVo® MULTIflex® Connections
20 Watts of Power
1-Year Warranty

KaVo® and MULTIflex® are registered trademarks of Kaltenbach & Voigt GmbH, Germany.
NSK® is a registered trademark of Nakanishi, Inc.

Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2022 Brasseler USA. All rights reserved. Printed in the U.S.A.

BRASSELER USA®