BRASSELER USA®

BY YOUR SIDE

FORZA PRO

PREMIUM ELECTRIC ATTACHMENTS


The Forza Pro series features a scratch-resistant glass coating for improved aesthetics and durability, a heat-blocking back cap that helps minimize the risk of patient burns, and the most extended warranty in the industry at three years! The full lineup includes a 1:5 High Speed, a 1:1 Latch, and a new 1:1 FG attachment, ideal for margin refinement.



- COMPATIBLE WITH ALL E-TYPE ELECTRIC MOTORS
- SCRATCH-RESISTANT GLASS COATING
- 5-PORT WATER SPRAY
- HEAT-BLOCKING BACK CAP
- SMALL HEAD SIZE
- 3-YEAR WARRANTY

- Order online at **BrasselerUSA.com** or call **800.841.4522**.
  In Canada call **800.363.3838**.

FORZA PRO





Forza™ Pro 1:5

Order Number: 5029811U0

1:5 FG
Fiber Optic
Head Size: Ø9.5 x H 14.5 mm
3-Year Warranty



Forza™ Pro 1:1 Latch

Order Number: 5029812U0

1:1 Latch
Fiber Optic
Head Size: Ø9.5 x H 12.4 mm
3-Year Warranty



Forza™ Pro 1:1 FG

Order Number: 5029813U0

1:1 FG
Fiber Optic
Head Size: Ø9.5 x H 13.8 mm
2-Year Warranty

▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2022 Brasseler USA. All rights reserved. Printed in the U.S.A.

BRASSELER USA®

REV 0 **B-5688**-2022-10-25