K0429 Peri-Implantitis Treatment Kit featuring C-Series™ TR™ Carbides

Developed with Dennis Tarnow, DDS, & Paul S. Rosen, DMD, MS

	Product Name	US No.	Order No.	Head Size ($^1/_{10}$ mm)	Cutting Length (mm)	Blade	Shank Type	Usage
	C-Series™ TR1 Long Flame Trimming & Finishing Carbide	-	TR1	014	8.0	F 12 Blade Fine	FG / XL Surg 30mm	Initial debriment of dental implant surface
	C-Series™ TR1UF Long Flame Trimming & Finishing Carbide	-	TR1UF	014	8.0	UF 30 Blade Ultra-Fine	FG / XL Surg 30mm	Final debriment of implants
	C-Series™ TR2 Long Flame Trimming & Finishing Carbide	-	TR2	023	8.0	F 12 Blade Fine	FG / XL Surg 30mm	Intial debriment of implant surface of crown contour makes it difficult for access to implant surface
	C-Series™ TR2UF Long Flame Trimming & Finishing Carbide	-	TR2UF	023	8.0	UF 30 Blade Ultra-Fine	FG / XL Surg 30mm	Final debriment of implant surface
	C-Series™ TR3 Football Trimming & Finishing Carbide	7404	TR2UF	014	3.1	F 12 Blade Fine	FG / XL Surg 30mm	Initial debriment of dental implant surface; Also provides osteoplasty to flatten bony contours where grafting isn't possible
	C-Series™ TR3UF Football Trimming & Finishing Carbide	7404	TR3	014	3.1	UF 30 Blade Ultra-Fine	FG / XL Surg 30mm	Final debriment of implant surface, bone or crown margin
	C-Series™ TR4 Football Trimming & Finishing Carbide	7408	TR3UF	023	4.1	F 12 Blade Fine	FG / XL Surg 30mm	Initial debriment of dental implant surface; Also provides osteoplasty to flatten bony contours where grafting isn't possible
	C-Series™ TR4UF Football Trimming & Finishing Carbide	7408	TR4	023	4.1	UF 30 Blade Ultra-Fine	FG / XL Surg 30mm	Final debriment of implant surface
	C-Series™ TR5 Flame Trimming & Finishing Carbide	-	TR5	016	3.5	F 12 Blade Fine	FG / XL Surg 30mm	Initial debriment of dental implant surface; Implantoplasty in deeper narrow defects
	C-Series™ TR5UF Flame Trimming & Finishing Carbide	-	TR5UF	016	3.5	UF 30 Blade Ultra-Fine	FG / XL Surg 30mm	Final debriment in base of defect for narrow peri-implant lesions

Recommended Speed: 40,000 rpm

▶ Visit our website at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.


BRASSELER USA®

©2023 Brasseler USA. All rights reserved.

K0429 Peri-Implantitis Treatment Kit featuring C-Series™ TR™ Carbides

Developed with Dennis Tarnow, DDS, & Paul S. Rosen, DMD, MS





Order No.	TR1	TR1UF	TR2	TR2UF	TR3	TR3UF	TR4	TR4UF	TR5	TR5UF
Head Size (1/10 mm)	014	014	023	023	014	014	023	023	016	016
Cutting Length (mm)	8.0	8.0	8.0	8.0	3.1	3.1	4.1	4.1	3.5	3.5
Friction Grip										
30 - XL FG Surg (30mm)	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg	XL FG Surg

Visit our website at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

BRASSELER USA®

©2023 Brasseler USA. All rights reserved.