K0433 Direct/Indirect Bonded Restorations System Usage Guide

CER BRD

	Product Name	US No.	Order No.	Head Size ($^1/_{10}$ mm)	Cutting Length (mm)	Grit/ Blade	Shank Type	Usage
	Round Diamond	-	801.31.023	023		M	31 - FG	Large onlay/overlay prep and large caries removal
	Round Diamond	-	801.31.018	018		M	31 - FG	Medium inlay/onlay/direct prep and medium caries removal
	Round Diamond	-	6801L.31.016	016		C	31 - FG	Preparation and caries removal with difficult access (e.g. deep margins, interproximal of neighboring tooth etc.)
	Round Diamond	-	801.31.014	014		M	31 - FG	Small direct prep and caries removal
	Round Diamond	-	801.31.010	010		M	31 - FG	Ultraconservative direct prep and caries removal
	Football Diamond	-	368.31.023	023	5.0	M	31 - FG	Lingual contouring (anterior) and occlusal reduction (posterior)
	Round End Taper Diamond	-	850.31.023	023	10.0	M	31 - FG	Veneer/crown/inlay prep with access
	Round End Taper Diamond	-	850.31.016	016	10.0	M	31 - FG	Veneer/crown/inlay prep (interdental)
	Round End Taper Diamond	-	850.31.014	014	10.0	M	31 - FG	Interdental prep/separation
	Round End Taper Diamond	-	6856.31.027	027	10.0	C	31 - FG	Large inlay/onlay/crown prep
	SabreCut Round End Cross-Cut Fissure	1558SC	1558SC	012	4.2	6-Blade	31 - FG	Amalgam/alloy removal
	Football Trimming & Finishing Carbide	7408	H379.31.023	023	4.2	F 12 Blade	31 - FG	Lingual (anterior) and occlusal adjustments (posterior) - Low speed
	Flame Trimming & Finishing Carbide	7901	H246.31.009	009	3.6	F 12 Blade	31 - FG	Fissure sculplting
	Flame Diamond	-	8274.31.016	023	3.5	F	31 - FG	Prep bevel and restoration margin finishing - Low speed
	Football Diamond	-	8368.31.023	014	5.0	F	31 - FG	Lingual (anterior) and occlusal prep finishing (posterior) and bevel
	Round End Taper Diamond	-	8856L.31.020	023	9.0	F	31 - FG	Veneer/crown/inlay prep margin finishing with access
	Round End Taper Diamond	-	8850.31.016	023	10.0	F	31 - FG	Veneer/crown/inlay prep margin finishing (interdental)
	Round End Taper Diamond	-	8850.31.014	014	10.0	F	31 - FG	Tight prep finishing (interdental)
	Round Diamond	-	8801.31.023	014		F	31 - FG	Large onlay/overlay prep margin finishing / Hollow chamfer with access
	Snap-On Mandrel		310.21				21 - RA Latch	Mandrel for flexible polishing/finishing discs
	VisonFlex™ Diamond Disc		934.21.180				21 - RA Latch	Interdental contact stripping (veneer preparation)





HANDPIECE SUGGESTIONS:

- Forza ELXT Electric Motor
- Forza M5 High Speed Electric Handpiece Attachment
- Forza M1 Low Speed Electric Handpiece Attachment





▶ Visit our website at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.


BRASSELER USA®

©2024 Brasseler USA. All rights reserved.

K0434 Composite/Ceramic Finishing System Usage Guide

CER BRD

	Product Name	Order No.	Head Size ($^1/_{10}$ mm)	Grit	Shank Type	Usage
	DiaLite™ Feather Lite™ Multi-Use Premium 2-Step Polisher	FLW14MZR.RA	140	Medium	21 - RA Latch	Prepolishing of ceramic/porcelain (step 1)
	Diacomp™ Point Multi-Use Premium 2-Step Polisher	DC1M21	040	Medium	21 - RA Latch	Prepolishing of ceramic/porcelain (step 2) and polishing of composite resin (step 1)
	Diacomp™ Feather Lite™ Multi-Use Premium 2-Step Polisher	FLDC5M	140	Medium	21 - RA Latch	Polishing of ceramic/porcelain (step 3) and composite resin (step 2)
	Diacomp™ Point Multi-Use Premium 2-Step Polisher	DC121	040	Fine	21 - RA Latch	Polishing of ceramic/porcelain (step 4) and polishing of composite resin (step 3)
	Diacomp™ Feather Lite™ Multi-Use Premium 2-Step Polisher	FLDC5	140	Fine	21 - RA Latch	Polishing of ceramic/porcelain (step 5) and composite resin (step 4)
	ComposiPro™ Brush	D209001.21	075		21 - RA Latch	Finishing of ceramic/porcelain (step 6) and composite resin (step 5), preferrably with water



HANDPIECE SUGGESTIONS:

- Forza ELXT Electric Motor
- Forza M5 High Speed Electric Handpiece Attachment
- Forza M1 Low Speed Electric Handpiece Attachment

▶ Visit our website at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

BRASSELER USA®

©2024 Brasseler USA. All rights reserved.