BRASSELER USA®


BRIOSHINE™ FEATHER LITE™ MP3™

MULTI-PURPOSE PURPLE POLISHERS

Embossed heads and color-coded for ease of identification

Diamond-impregnated polyurethane rubber with a unique shape and flexibility; provides access to areas other shapes can't

Conveniently sterile for single-patient use; ready to use, disposable



Starts as a knife-edge used as a finishing disc then wears into starburst-shaped ribs; adapts to all tooth surfaces

Gold plastic latch mandrel for use with low speed handpieces*

An Innovative Polisher That Hits All The Right Notes

MP3 is multi-purpose and can be used on **all restorative materials** to become your go-to polisher. Indications include:

- Orthodontic removal of residual aligner attachments
- Removal of excess cement
- Smooths enamel scratches
- Removal of surface stains

Additional indications listed at
BrasselerUSADental.com/MP3



Order No.	MP3C	MP3M	MP3F
Head Size (1/10 mm)	100	100	100
Grit	Coarse	Medium	Fine
Color	Dark Purple	Light Purple	Gray
21 - RA - Right-Angle Latch (22 mm)			
Qty/Pack	25	25	25

Recommended speed: 5,000-8,000 rpm
Final Polish: 2,000-4,000 rpm

**Brasseler Low Speed Handpiece Recommendations:*
- *Forza Pro 1:1 Latch Electric Attachment*
- *Forza LS Latch Air Attachment*

▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2023 Brasseler USA. All rights reserved. Printed in the U.S.A.