BRASSELER USA®

CORE BUILD-UP MATERIAL

DUAL-CURE WITH FLUORIDE

When there is insufficient tooth structure available for a fixed prosthesis, Brasseler USA's Core Build-Up Material provides the added strength and stability needed for long-term clinical success.



CARTRIDGES

1	Cartridge
25	Mixing Tips
25	Intra-Oral Tips
1	Mixing Pad

Shade	Size	Order No.
Enamel	50g	5031045U0
White	50g	5031047U0



Mixing Tips


20/Pack 5031049U0

Intra-Oral Tips


20/Pack 5031297U0

SYRINGES

2	Syringe
20	Mixing Tips
20	Intra-Oral Tips



Shade	Size	Order No.
Enamel	10g	5031044U0
White	10g	5031046U0



Mixing Tips

20/Pack 5031048U0

Intra-Oral Tips

20/Pack 5031296U0

- **DUAL-CURE**
  Provides confidence of a self-cure with option of light-curing to reduce chair time
- **CUTS LIKE DENTIN**
  For easy shaping
- **RADIOPAQUE**
  Easily distinguishable in radiographs
- **CONTAINS FLUORIDE**
  Helps prevent secondary decay



- **STACKABLE VISCOSITY**
  Allows build-up without a matrix or core form

Indications:

- *Cementation of posts*
- *Composite for crown build-ups*
- *Repairing and preventing split roots*
- *Integral post and core procedures in one visit*
- *Enamel shade may also be used as a posterior restorative material*

POST & CORE SOLUTIONS

BRASSELER USA®

ENDOSEQUENCE® POSTS *(5/Pack)*



White and translucent for highly esthetic restorations and extremely strong with modulus of elasticity like that of dentin, which results in less stress on the tooth.

Tip ID	Tip Size	.04 Taper	.06 Taper
Yellow	.50	5013508U0	5013512U0
Red	.70	5013509U0	5013513U0
Blue	.90	5013510U0	5013514U0
Green	1.0	–	5013515U0
Green	1.1	5013511U0	–

ENDOSEQUENCE® BC LINER™ *(5 mL)*

Provides a protective barrier between the dental pulp and restorative material.


BC Liner™ 5027169U0



BC Liner™ Blue 5027706U0

CROWN PREPARATION DIAMONDS

Distinctive high-quality diamonds and handpieces to make this step in the procedure effortless and precise every time.



US No.	–	–	KS1	–	–
Order No.	**5856**	**6856DC**	**35005**	**WS8M1.0**	**5847KR**
Head Size (1/10 mm)	016	016	012	018	016
Cutting Length (mm)	8.0	8.0	7.0	8.0	8.0
Grit	Super-Coarse	DuraCut®	Super-Coarse	Medium	Super-Coarse
Shank Type	FG, FG/SS	FG	FG, FG/SS	FG	FG, FG/SS

FORZA™ PRO HIGH SPEED ELECTRIC ATTACHMENT

1:5 Speed Increasing
Fiber Optic
Head size: 14.5mm x 9.5mm
3-Year Warranty



Order No. 5029811U0

MAX ETCH 35% PHOSPHORIC ACID

Gel-like viscosity ideal for the selective-etch technique, helping to prevent sensitivity and to protect gingival tissue.

5 mL 345503-B



Not Available in Canada

MPA UNIVERSAL



Provides consistent and high bond strengths regardless of placement technique.

MPa Universal	Size	Qty.	Order No.
Light-Cure Adhesive	0.1 mL	50	267000-B
Light-Cure Adhesive	5 mL	1	267001-B
Dual-Cure Catalyst	5 mL	1	267002-B

Not Available in Canada

▶ Order online at **BrasselerUSA.com** or call **800.841.4522**.
In Canada call **800.363.3838**.

©2024 Brasseler USA. All rights reserved. Printed in the U.S.