BRASSELER USA®



K0444 Apa: Aesthetic System

In collaboration with Michael Apa, D.D.S.

	Product Name	Order No.	Head Size (1/10 mm)	Cutting Length (mm)	Shank Type	Grit	Recommended Speed (RPM)	Application	Order Quantity
	Big Daddy	5856L	025	9.0	31-FG (19 mm)	Super-Coarse	200,000	· Bulk reduction for veneer preparation.	
	The Chamfer	6850	016	10.0	31-FG (19 mm)	Coarse	200,000	· Refinement of margins for veneer preparation.	
	Green and Gold	2847KR	018	8.0	31-FG (19 mm)	CRF™	200,000	· Bulk reduction for crown preparation.	
	Sculptor	859	010	11.0	31-FG (19 mm)	Medium	200,000	· Adjusting facial of temps.	
	Flame	889	009	3.5	31-FG (19 mm)	Medium	200,000	· Interproximal reduction for temp removal. · Beveling/refinement between papilla.	
	Football	8379	021	4.2	31-FG (19 mm)	Fine	200,000	· Occlusal reduction for bite check.	
	IPS e.max® Football	6379	023	4.2	31-FG (19 mm)	Coarse	200,000	· Bulk reduction for bite check.	
	Mosquito	8392	016	5.0	31-FG (19 mm)	Fine	200,000	· Refinement of margins for veneer insert. · Refinement of interproximal for temps.	
	T-Polish	859UF	016	10.0	31-FG (19 mm)	Ultra-Fine	200,000	· Finishing interproximal provisionals.	
	No Band	847KR	016	8.0	31-FG (19 mm)	Medium	200,000	· Shoulder beveling for crown preparation.	
	Carbide	1558 (H031RS)	012	4.2	31-FG (19 mm)	–	200,000	· Alloy removal.	
	Inlay Bur	6845KR	018	4.0	31-FG (19 mm)	Coarse	200,000	· Inlay preparation.	
	Corgi	845KR	025	4.0	31-FG (19 mm)	Medium	200,000	· Large/fat inlay preparation.	
	Mandrel	310.21	–	–	21-RA (22 mm)	–	–	–	
	Spiral*	T5856	016	8.0	31-FG (19 mm)	Super-Coarse	200,000	· Removal of zirconia.	

**Not included in K0444 kit*

IPS e.max® is a registered trademark of Ivoclar Vivadent®

Customer Number

Address

Business Name or Practice Name

Practitioner Name

City State Zip

Email Address

( )

Business Phone

To order, visit **Shop.BrasselerUSA.com** or call **800.841.4522**.
In Canada, call **800.363.3838**.

BRASSELER USA®

©2026 Brasseler USA. All rights reserved.

26-1561 | 2026-07-17 REV0